Handleiding glazen pH-elektrode

Prepareren

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Instructions Glass pH Electrode

Preparation

Unscrew the electrode from the "storage bottle" which contains 50/50 mixture of pH 4 buffer and 4M KCl. Carefully slide the bottle and the O ring from the electrode. Rinse the glass ball and the ceramic part with tap water. If there are air bubbles in the glass ball, shake the electrode until the ball is completely filled with fluid. **Keep the bottle with fluid and the O ring for storage.**

Calibration

In general: see your system machinery for the instructions for use.

1. Rinse the tip of the electrode with distilled water and stick it in a gauging buffer with pH 7.0. Stir the electrode in the buffer fluid for better contact and to remove any air bubbles.
2. Give the electrode (a few minutes) time to stabilise. Adjust the pH meter with the gauging button to "7.0".
3. Remove the electrode from the buffer fluid, rinse well with distilled water, shake off the excess water and stick the electrode in buffer fluid with pH 4.0.
4. Again, give the electrode time to stabilise and adjust the pH meter to the right pH using the slope button. **The electrode should always be rinsed in between with distilled water to prevent dilution.**
5. For very accurate measurements, it is advisable to repeat steps 1 through 4 to ensure complete stabilisation of the electrode.

Storage

1. For short-term storage it is advisable to dip the electrode in the 50/50 mixture of pH 4 buffer and 4M KCl. **Do not store in distilled water!**
2. For long-term storage, the plastic cap or storage bottle should be filled with 50/50 mixture of pH 4 buffer and 4M KCl. Stick the electrode in the fluid. Remove any excess fluid. **Do not store in distilled water!**

Cleaning

The functioning of used electrodes can be improved by one of the following cleaning procedures:

1. General: It is advisable, to regularly clean the glass ball and the ceramic part of the electrode with a solution of fluid detergent in warm water. Use a soft toothbrush or a tissue, dipped in the soap solution. Rinse the tip well with distilled water. Then dip the electrode in 50/50 mixture of pH 4 buffer and 4M KCl for at least 30 minutes before using it again. **Do not store in distilled water!**
2. Anorganic deposit: Dissolve the deposit by for a few minutes dipping the electrode in diluted 1:1 nitric acid or phosphoric acid 1:4; Nitric: add 1 part acid to 1 part water, Phosphoric: add 1 part acid to 4 parts water. Then rinse the ball and ceramic part well with tap water. Continue with point 1, as described above. **Add the acid to the water.**

Remarks

1. Do not let the electrode dry out.
2. Torn or broken electrodes cannot be repaired, replace the electrode immediately.
3. Inspect the cable and plug to ensure that the insulation is still intact and that there are no signs of corrosion or contamination of the metal parts. Broken cable insulation and snapped or knotted cables can give rise to bad functioning of the electrode.

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