



## SAFETY DATA SHEET

### 1. Product and Company Identification

Product Name	Electrode Cap Storage Solution
Product Number	R001000
Product Use	Calibration Solution
Manufacturer	Van London Company
Address	10540 Rockley Road, Houston, Texas 77099
Telephone	832-456-6641
Emergency Phone	832-456-6641

### 2. Hazards Identification

#### Emergency Overview

##### OSHA Hazards

Acute aquatic toxicity (Category 3), H402  
Chronic aquatic toxicity (Category 3), H412

##### HMIS Classification

**Health Hazard:** 1  
**Flammability:** 0  
**Physical Hazards:** 0

##### NFPA Rating

**Health Hazard:** 1  
**Fire:** 0  
**Reactivity Hazard:** 0

##### Potential Health Effects

**Inhalation:** No effects anticipated  
**Skin:** No effects anticipated  
**Eyes:** No effects anticipated  
**Ingestion:** No effects anticipated

### 3. Composition/Information on Ingredients

Component	CAS Number	Weight %
Potassium hydrogen phthalate	877-24-7	0.5
Potassium chloride	7447-40-7	15
Deionized water	7732-18-5	84.5

### 4. First Aid Measures

<b>Eye</b>	Flush eyes with water as a precaution.
<b>Skin</b>	Wash off with soap and plenty of water.
<b>Inhalation</b>	Move person into fresh air. Not breathing, give artificial respiration.
<b>Ingestion</b>	Give large quantities of water. If you feel unwell, contact a physician.

## 5. Fire Fighting Measures

### Suitable

**extinguishing media** Water, CO<sub>2</sub>, Dry chemical or foam

### Hazardous

**combustion products** Hydrogen chloride gas, Potassium oxides.

### Special

### protective

### equipment for

### firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

## 6. Accidental Release Measures

### Personal precautions

Avoid contact with eyes. Ensure adequate ventilation.

### Environmental precautions

Do not let product enter drains. Discharge into the environment must be avoided.

### Methods and materials for containment and cleaning up

Soak up with inert absorbent material and dispose of as hazardous waste. Keep in suitable, closed containers for disposal.

## 7. Handling and Storage

### Handling

Avoid contact with eyes. Wash thoroughly after handling. Avoid inhalation of vapour or mist.

### Storage

Keep container tightly closed and upright in a dry and well-ventilated place. Store at ambient or lower temperature. Protect against physical damage.

## 8. Exposure Controls and Personal Protection

### Exposure Limits

Contains no substances with occupational exposure limit values

### Engineering Controls

Maintain general industrial hygiene practices when using this product.

### Personal Protective Equipment

#### Respiratory Protection

Respiratory protection not required. For nuisance exposures use type OV/AG (US) or type ABEK (EU EN 14387) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

#### Hand Protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Wash and dry hands after handling.

#### Eye Protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).

#### Skin and body protection

Impervious clothing.

#### Hygiene measures

General industrial hygiene practice.

## 9. Physical and Chemical Properties

### Appearance

Clear, colorless liquid

### Odor

None

<b>Odor Threshold</b>	No data available
<b>pH</b>	4.01
<b>Freezing Point</b>	0°C
<b>Boiling Point</b>	100°C
<b>Flash Point</b>	Not flammable
<b>Evaporation rate</b>	No data available
<b>Flammability</b>	No data available
<b>Lower explosion limit</b>	No data available
<b>Upper explosion limit</b>	No data available
<b>Vapor pressure</b>	No data available
<b>Solubility</b>	Soluble in water
<b>Vapor density</b>	No data available
<b>Relative density</b>	No data available
<b>Partition coefficient:</b>	
<b>n-octanol/water</b>	No data available
<b>Auto-Ignition temperature</b>	No data available
<b>Decomposition temperature</b>	No data available
<b>Density</b>	No data available
<b>Viscosity</b>	No data available

#### 10. Stability and Reactivity

<b>Chemical Stability</b>	Stable under recommended storage conditions.
<b>Possibility of Hazardous reactions</b>	No data available
<b>Conditions to avoid</b>	Extreme temperatures
<b>Materials to avoid</b>	Strong oxidizing agents
<b>Hazardous decomposition products</b>	None

#### 11. Toxicological Information

<b>Signs and Symptoms of Overexposure</b>	To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.	
<b>Acute Effects</b>	<b>Eye Contact</b>	No data available
	<b>Skin Contact</b>	No data available
	<b>Inhalation</b>	No data available
	<b>Ingestion</b>	No data available
	<b>Germ Cell Mutagenicity</b>	No data available
<b>Target Organ Effects</b>	No data available	
<b>Chronic Effects</b>	No data available	
<b>Medical Conditions Aggravated by Exposure</b>	No data available	
<b>Acute Toxicity Values</b>	<b>Oral LD50</b>	None reported
	<b>Inhalation LC50</b>	None reported
	<b>Dermal LD50</b>	None reported
<b>Carcinogenicity</b>	IARC	No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

**Reproductive Toxicity** No data available

**Teratogenicity** No data available

**Specific target organ toxicity - single exposure** - No data available

**Specific target organ toxicity - repeated exposure** - No data available

**Aspiration Hazard** No data available

**Potential Health Effects**

<b>Inhalation</b>	No data available
<b>Ingestion</b>	No data available
<b>Skin</b>	No data available
<b>Eyes</b>	No data available

**12. Ecological Information**

**Toxicity** No data available

**Persistence & degradability** No data available

**Bioaccumulative potential** No data available

**Mobility in soil** No data available

**Other adverse effects** No data available

**13. Disposal Considerations**

**Disposal method** Offer surplus and non-recyclable solutions to a licensed disposal company.

**Contaminated packaging** Dispose of as unused product.

**14. Transport Information**

**DOT (US)** Not dangerous goods

**IMDG** Not dangerous goods

**IATA** Not dangerous goods

**15. Regulatory Information**

**OSHA Hazards** No known OSHA hazards

**SARA 302 Components** No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

**SARA 313 Components** This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

**SARA 311/312 Hazards** Chronic Health Hazard

**Massachusetts Right To Know Components** No components are subject to the Massachusetts Right to Know Act.

**Pennsylvania Right To Know Components**

	CAS-No.	Revision Date
Potassium chloride	7447-40-7	
Water	7732-18-5	

**New Jersey Right To Know Components**

	CAS-No.	Revision Date
--	---------	---------------

Water	7732-18-5
Potassium hydrogen phthalate	877-24-7
Potassium chloride	7447-40-7

**California Prop. 65 Components**

This product does not contain any chemicals known to the State of California to cause cancer, birth defects, or any other reproductive harm.

**16. Other Information**

**Issuing Date** 12/16/2015

**Revision No.** 1.1

**Disclaimer** The information contained herein is accurate to the best of our knowledge, but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. Van London pHOenix Company makes no warranty of any kind, expressed, or implied, concerning the safe use of this material in your process or in combination with other substances. Van London pHOenix Company shall not be held liable for any damage resulting from handling or from contact with the above product.