



VAN LONDON
Est. 1961



"When Accuracy Matters"

When Van London Company was founded in 1961, Gregory S. Van London had a special dedication to providing products and services which would either meet or exceeded customers' expectations. Today, we continue his philosophy with conviction. Consumer satisfaction is our greatest ambition, and 30-plus years of manufacturing and applications experience assures it. Working with our customers creates a never-ending process of improvement and innovation for all aspects of our business.

pH · ORP · ISE · DO · Conductivity · Custom · Direct Replacements

Why Choose VL-PC ?

- Quality Product
- Customized Service
- Stability in the Market Place
- Premium Customer Service
- Dedicated R&D
- Reliability
- Experienced Tech. Support
- All Associated Accessories
- Fast Turn-around Service
- Competitive Pricing
- Custom Manufacturing
- Private Label
- Flexible Stocking Options
- Direct Replacement For Most Major Brands
- Extensive quality control
- Technical Application Assistance




· Full Machine Shop
State of the art equipment



· Newly Refurbished 40,000 ft2 (3 716 m2)
Manufacturing Plant
· Bar Coded Inventory

Developers
of the
Twist-Lock™



- New pH Glass formula
- 0-14 pH
- 10 Seconds(95% Response)
from pH 7.01 to pH 4.01
- Glass Membrane Impedance
20 Mohm at 25°C



- Tuff-Tip ®
- Rugged pH Tip
- Easy to Clean
- Fast Response
- Non-Clogging Reference Junction
- Seal Cross Linked Polymer Reference



Applications



Industrial



Bio Tech



Laboratory

Van London-pHOenix Co.
6103 Glenmont Drive
Houston, Texas 77081
Local: (713) 772-6641
Toll Free: (800) 522-7920
Web: www.VL-PC.com
e-mail: info@VL-PC.com
We Speak Spanish & Vietnamese



Made in the USA

U.S. Small Business Administration