



Introducing our new pH glass formula with unparalleled performance. This new pH glass formula,GXV, can replace all other formulas, streamlining inventories and eliminating wasted time trying to select the best electrode for the job. The GXV formulation is perfect in low ionic strength solutions as well as high pH applications and is specially formulated to have lightening fast response time with superior accuracy and repeatability. The sample pH is registered in less than 10 seconds, the result is correct to within 0.02 pH units, and the pH value is reproducible over the entire pH range from 0 to 14 pH.

For the best pH glass formula on the market, Just Think Green!



6103 Glenmont Drive Houston, Texas 77081 Local: (713) 772-6641 Toll Free: (800) 522-7920 Web: www.VL-PC.com e-mail: info@vl-pc.com Se Habla Español

Van London-pHoenix Co.