CALIBRATION & OTHER SOLUTIONS

1/11



6103 GLENMONT • HOUSTON, TEXAS • 77081 800/522-7920 • FAX 713/772-4671 • 713/772-6641 www.VL-PC.com

ORDERING INFORMATION

To build your own part numbers, replace X's with selections listed below.

921XXXX	: Calibration & Reference Solutions
0012	475ml of 4 pH buffer solution, NIST-traceable
0013	475ml of 7 pH buffer solution, NIST-traceable
0014	475ml of 10 pH buffer solution, NIST-traceable
0016	•
0230	475ml of 4 Molar KCI (gelled RFH)
0082	,
0102	
0008	475ml of Cap Solution (pH electrode storage solution)
0146	475ml of Conductivity Standard 146.9 microSiemens
0717	
1412	475ml of Conductivity Standard 1,412 microSiemens
2765	475ml of Conductivity Standard 2,765 microSiemens
6667	
0800	475ml of Anti-Fouling Solution
2011300	475ml of D.I. water (rinse solution)
2011300-1G	1 gallon size of D.I. water (rinse solution)
9210012-1G	1 gallon size of 4 pH buffer solution, NIST-traceable
9210013-1G	1 gallon size of 7 pH buffer solution, NIST-traceable
9210014-1G	1 gallon size of 10 pH buffer solution, NIST-traceable

- 4, 7, and 10 pH buffer solutions are color coded for visual convenience.
- Special solutions and Ion Selective solutions are available. Consult factory for options.

FOR ADDITIONAL PRODUCT INFORMATION, CALL 713/772-6641 OR 800/522-7920 E-mail: info@VL-PC.com www.VL-PC.com

[Buffer.indd]