

WIRING ADAPTERS

2/07



6103 GLENMONT • HOUSTON, TEXAS • 77081
800/522-7920 • FAX 713/772-4671 • 713/772-6641
www.vl-pc.com

ORDERING INFORMATION

- To build your own part numbers, replace **X**'s with selections listed in the drop-down menu.

2010XXX - XX	: VL-PC Common Wiring Adapters
	Lead Length
	Specify lead length in feet.
	Adapter Styles
400	SL to SL (2 conductor)
520	SL to SL (coax)
590	SL/SL to SL/SL (4 conductor with drain) †
690	SL to SL (6 conductor with drain)
610	SL to SL (coax with 2 conductor)
507	SL to Hanked (4 conductor with drain) †
531	SL to Hanked (coax with 2 conductor)
691	SL to Hanked (6 conductor with drain)
480	FBNC to MBNC (coax)
430	FBNC to SL (coax)
506	FBNC/Hanked to SL (coax with 2 conductor)
530	FBNC/SL to MBNC/SL (coax with 2 conductor)
650	FBNC/SL to SL/SL (coax with 2 conductor)
680	FBNC/SL to MBNC/Hanked (coax with 2 conductor)
655	FBNC/SL to Hanked (coax with 3 conductor) **
512	SL to MBNC (coax)
670	SL/SL to MBNC/SL (coax with 2 conductor)
620	FBNC/Male Molex to Q/6 pin (for L&N)
681	FBNC/Male Molex to MBNC/Hanked (coax with 2 conductor)
710	FBNC/Male Molex to SL (coax with 2 conductor)
730	FBNC/Male Molex to Hanked (coax with 2 conductor)
350	MBNC to Alligator Clips (coax)
456	F4-pin Military to Hanked
452	F4-pin Military to male BNC/Hanked
662	SN6 to "Specify" (3mm Coax)
VP0	Female VP to "Specify"
VP1	Female VP with Solution Ground to "Specify"


- † Standard cable for use with Van London-Phoenix Co. Conductivity Sensors.
- * Consult the factory for other available options.
- * All coax is constructed of special (graphite) cable.
- ** For use with Solution Ground sensors

CONNECTORS

ORDERING INFORMATION

- To build your own part numbers, first determine the connector style required and then replace X's with selections from within the parentheses listed below. (ie. . .79B5760-10X)

male BNC "B"



female BNC with male Molex "WM"




2-Pin female Molex "T"



male BNC with female Molex "BT"




female BNC with Spade Lug "WS"



Standard Plug "S"




male BNC with Spade Lug "BS"



female BNC with Hanked "WH"




4-Pin Military (on electrode) "F"




male BNC with Hanked "BH"




Spade Lug "SL"




Ingold (Mettler-Toledo) Cap (on electrode) "S8"




male BNC with Miniature Audioplug "BMA"




Hanked "H"




Detachable Lead (on electrode) "DL"




female BNC "W"



2-Pin male Molex "M"



Ingold (Mettler-Toledo) Cap (on electrode) "S7"



**FOR ADDITIONAL PRODUCT INFORMATION,
CALL 713/772-6641 OR 800/522-7920
E-mail: sales@VL-PC.com
www.VL-PC.com**

* Consult factory for other styles available.

[WiringAdapters.pmd]